

Fixed Volume and Triple Volume Micro Pipette

Nichipet F&V

- The light and slender design along with the smooth movement of the plunger reduces the operator's fatigue.
- Built-in tip ejector.
- PTFE (Fluorine resin) is used in the airtight chamber of the instrument, so the instrument keeps its air tight seal and precise reproducibility for long hours.



Nichipet F

- ▷ Suitable when only one volume is needed.
- ▷ A wide range of standard volumes from 10µL – 1000µL.
- ▷ Other volumes within the range available for special order.

Nichipet V

- ▷ Three different volume models.
- ▷ Volume change by simply turning the knob.

Specifications

Type	Cat. No	Volume setting	Volume range(µL)	Accuracy(%)	Precision(%)	Usable tips
F	00-NPF-10	Fixed	10	±1.2	≤0.4	BMT2-SS/SSR, FSS/FSSR
	00-NPF-20		20	±0.8	≤0.3	BMT2-SG/SGR, SE/SER, FSG/FSGR
	00-NPF-25		25	±0.8	≤0.3	BMT2-SG/SGR, SE/SER, FSG/FSGR*, FSE/FSE
	00-NPF-50		50	±0.6	≤0.3	
	00-NPF-100		100	±0.6	≤0.3	
	00-NPF-200		200	±0.6	≤0.2	BMT2-SG/SGR, SE/SER, FSE/FSE*, FAG/FAGR
	00-NPF-500		500	±0.5	≤0.2	BMT2-LG/LGR, LE/LE
	00-NPF-1000		1000	±0.5	≤0.2	BMT2-LG/LGR, LE/LE, FLG/FLGR, FLE/FLER
V	00-NPV-S	Three variable settings	10/20/50	±1.2 – ±0.7	≤0.6 – ≤0.3	BMT2-SG/SGR, SE/SER, FSG/FSGR
	00-NPV-M		50/100/200	±0.7 – ±0.5	≤0.4 – ≤0.3	BMT2-SG/SGR, SE/SER, FSE/FSE*, FAG/FAGR
	00-NPV-L		200/500/1000	±0.7 – ±0.5	≤0.3 – ≤0.2	BMT2-LG/LGR, LE/LE, FLG/FLGR, FLE/FLER

*The AC and CV values are the values by use of the disposable tips described in the catalogue.
 *See P15-16 for tip selection.

Accessories

Stand



Cat. No 00-MLT-STD2(ABS)
 Storing capacity 6pcs
 Compatibility Nichipet Premium
 Nichipet EXII
 Nichipet EX PlusII
 Nichiryo Le



Cat. No 00-BM-STD(Acrylic resin)
 Storing capacity 6pcs
 Compatibility Nichipet Premium
 Nichipet EXII
 Nichipet EX PlusII
 Nichipet F&V




Cat. No 00-NLE-STD(ABS)
 Storing capacity 1pcs
 Compatibility Nichipet Premium
 Nichipet EXII
 Nichipet EX PlusII
 Nichiryo Le



Cat. No 00-7000S(PVC)
 Storing capacity 2pcs
 Compatibility Nichipet 7000

Tray



Cat. No 00-88-TY
 Q'ty 50pcs
 Compatibility Nichipet 7000

Filter

For NPP-1000

Cat. No ... 00-NP-0601000
 Q'ty 10pcs

For NPP-5000,10000

Cat. No ... 00-NP-0615000
 Q'ty 100pcs

For NPXII-1000

Cat. No ... 00-NX2-0601000
 Q'ty 10pcs

Tube Opener



Cat. No 00-NTO-001
 Q'ty 3pcs
 Compatibility Nichipet Premium
 Nichipet EX II
 Nichipet EX Plus II

For NPLOII-1000

Cat. No ... 00-LO2-0601000
 Q'ty 10pcs

For NPXII, NPLOII-5000,10000

Cat. No ... 00-NX2-0615000
 Q'ty 100pcs