

# EO GAS STERILIZER

2010  
NEW PRODUCT



## FEATURES

- 1) It is very easy to check all of working situation by Use-friendly controller
- 2) When you use the gas, it counts automatically, so that the gas level is convenient to check.
- 3) Safety is guaranteed by 12 self-checking system.
- 4) The best sterilization condition is possible by sprayed gas.
- 5) The equipped pump has long life time and strong for sterilization.

Program	Pressure	Temperature	Sterilizer	Aeration
1	1 Kg/cm <sup>2</sup>	38°C	5 Hour	12 Hour
2	1 Kg/cm <sup>2</sup>	55°C	3 Hour	9 Hour
3	User Select			

## SPECIFICATION

Spec. / Model		K-MKE030	K-MKE050
Overall Size (W x D x H, mm)		535 x 700 x 430	619 x 712 x 520
Size	Chamber Size (mm)	Ø261.4 x 540(L)	Ø350 x 510(L)
	Chamber Capacity	30ℓ	50ℓ
	Chamber Material	Stainless steel	
Tray	Size	515 x 262.5(mm) x 1ea	480 x 360(mm) x 1ea
	Tray Material	Stainless steel	
Using Gas		Eo(20%), CO(80%)	
Vacuum Pump		Oil-less Diaphragm Pump (40L/Min)	
Pressure		1 Kg/cm <sup>2</sup>	
Volt		AC 220V 50/60Hz	
Power Supply		800W	900W
Controller System		8Bit Microprocessor	
Weight		55Kg	68Kg
Cart (W x D x H, mm)		610 x 730 x 800	
Accessory		Gsa Tank, Sealing Machine, Roll Bag.	