

# FUME HOOD



## FEATURES

- 1) Cup-board is separatable for convenience.
- 2) The valve for water supplying is equipped.
- 3) Both of power supply for 110V, 220V are available.
- 3) Various utility valves are equipped.
- 4) The tempered glass is used for resercher's safety.

## DISTINCTION

1. Knob



Valve type.

2. Connector



Various Connector.

3. Water Supplier



Adjustable.

4. Conent



Both of 100V & 220V.

5. Valves



Convenient Valves

6. Tempered Glass



Safe Tempered Glass.

## SPECIFICATION

Spec. / Model		K-FH120	K-FH150	K-FH180
Size	External (WxDxH, mm)	1200x800x2300	1500x800x2300	1800x800x2300
Utility Device		Air Cock, Gas Cock, Water Cock x 1ea		Air Cock, Gas Cock, Water Cock x 2ea
Work Top		110 / 220 V consents, Lamp x 1ea		
Fan		Ceramite & Stainless steel (optional)		
Material		Φ200 (Built-In)		
		SUS #304 Stainless Steel		
		Steel Powder Coating		
Power		220V 50/60 Hz		

\* Customized model is possible

\*\* Above speficiation can be modified without notice.