

# HOT PLATE STIRRER



## FEATURES


- 1) Heating and stirring can operate separately.
- 2) Easy controls of heating temperature and stirring speed.
- 3) The ceramic coated aluminum top plate possessing a strong chemical durability and a quick heating ability has an excellent durability for long time use.
- 4) The anticorrosive stainless body
- 5) Ultra strong magnetic force prevents the spin bar from jumping.
- 6) Heavy duty motor for stable long time use.

## SPECIFICATION

Spec. / Model	K-HYHS11	K-HYSH12	K-HYSH13
Type	Hot plate & Magnetic stirrer	Hot plate	Magnetic stirrer
Temp range	Up to 320 ° C	Up to 320 ° C	-
RPM	2000	-	2000
Watt	680	650	30
Plate (WxDxH,mm)	200x200	200x200	200x200
Dimension (WxDxH,mm)	210x220x120	210x220x120	210x220x120
Electric Requirements	220VAC. 50/60Hz	220VAC. 50/60Hz	220VAC. 50/60Hz

\* Above specification can be modified without notice.

# MULTI ROOM SHAKER INCUBATOR



Spec. / Model	K-HB201MS2R	K-HB201MS4R
Type	2 Room 2 Shaker	4 Room 4 Shaker
Dimension (IN)	650x600x500(H)mm x 2 Room	650x600x500(H)mm x 4 Room
(OUT)	830x870x1750(H)mm	1800x870x1750(H)mm
Holder Plate(Total)	250mlx50EA	250mlx100EA
Temp. Range	0°C to 60°C ±0.2°C	
Temp. Controller	Digital P.I.D Controller	
Speed Controller	Digital Display Feed back Speed Controller	
Timer	Digital Timer 99hr 59min. or Infinity	
Power	2.5KW	4.5KW
Spec. / Model	K-HB201MS2	K-HB201MS4
Shaker	2 Shaker (1 Room)	4 Shaker (1 Room)
Dimension (IN)	650x600x1100(H)mm	1300x600x1100(H)mm
(OUT)	830x870x1700(H)mm	1480x870x1700(H)mm
Holder Plate(Total)	250mlx50EA	250mlx100EA
Temp. Range	0°C to 60°C ±0.2°C	
Temp. Controller	Digital P.I.D Controller	
Speed Controller	Digital Display Feed back Speed Controller	
Timer	Digital Timer 99hr 59min. or Infinity	
Power	1.8KW	2.5KW

\* Above specification can be modified without notice.