

CIRCULATION WATER BATH (HEAD TOP)

FEATUREA

- 1) Low height and vibration provide stable working
- 2) The accurate control by Digital controller
- 3) Various Flask and Spring Wire Rack are provided.
- 4) Belt and pulley is not used so easy maintenance and longer life time is possible.



DISTINCTION

1. Spring Wire Rack.



Various flask holder and spring wire rack can be used on customer needs.

1. Low Height



The low height provide stable working.

1. Flask Holder



Many kinds of flask holder.

SPECIFICATION

Spec. / Model		K-OS300	K-OS400	K-OS500
Motion		Orbital		
Temp	Range	25 to 300 rpm(up to 350rpm : high speed type)		
	Accuracy	±1%		
	Increment	1 rpm		
Time	Range	59min 59sec or Continuous		
	Accuracy	±1%		
	Increment	1 minute		
Controller		Microprocessor Digital		
Display		4 Digit LED , 2 Display		
Motor		Plate type B/L DC Motor		
Drive system		Beltless Type		
Stroke		20mm		
Size	Plarform (W x D mm)	330x330	470x470	730x570
	External (W x D x H, mm)	353x400x160	493x540x160	753x640x160
Electric Requirements		220VAC 50/60Hz		

** Above speficiation can be modified without notice.