

DIGITAL OVEN



FEATURES

- 1) Specially designed low-noise and low-vibration fan motor. (Forced convection model only)
- 2) The grounded circuit prevents a short circuit and overflow of electricity.
- 3) Equipped overheat protection sensor.
- 4) Accurate temperature controls by the Digital PID controller
- 5) High temperature ceramic wool & glass wool optimize insulation and heat loss.
- 6) Threefold tempered safety glass windows for excellent insulation and easy observation
- 7) Set time/temperature function
- 8) Stable standing with uniquely designed leveling foot.
- 9) Equipped door damper for airtight closeness and smooth operation.
- 10) Well-ventilated double doors prevent exterior heating and provide the safety for the operators.
- 11) Double doors prevent exterior heating and provide the safety for the operator.

DISTINCTION

1. The wide durable threefold glass window



Threefold easy to view observation Window with excellent insulation.

2. Easy to use door handle.



Door open & close easy operation.

3. The shelves easy to move



Well ventilated height adjustable stainless steel shelves, easy to clean up.

4. Overheat and over current protection



For safety, overheat, over current Protection installed.

5. PID controller



Strong noise resistant, various functions, easy reset, wonderful after sales service.

6. Slip-resistant casters.



Uniquely designed height adjustable, slip-resistant casters.