


## SPECIFICATION

Spec./ Model	Natural Convection	K-DO42	K-DO81	K-DO150
	Force Convection	K-CO42	K-CO81	K-CO150
Type		1 door		
Convection Type		Natural convection		
Temp	Range	Max 250 ° C		
	Controller	Digital P.I.D controller		
Size	Internal (W x D x H, mm)	350x350x350	450x400x450	600x500x500
	External (W x D x H, mm)	520x525x720	620x575x820	770x675x870
Safety	Electrical	No fuse breaker		
		Over Temp limit.		
Material	Internal	Stainless steel		
	External	Steel, Powder coating		
	Shelves (W x D, mm)	SUS #304 Wire Electro-Polished		
		327x300, 2ea	427x350, 2ea	577x450, 3ea
	Heater	Sheath heater		
	Insulation	Ceramic wool+ Glass wool		
	Door Packing	High Temp. grade foamed silicone rubber		
Windows	Tempered safety glass, Three fold			
Electric Requirements		220VAC. 50/60Hz		

\* Customized model is possible.

\*\* Above specification can be modified without notice.

## AIR CIRCULATION



**FORCED CONVECTION OVEN**

**K-CO Series**



**NATURAL CONVECTION OVEN**

**K-DO Series**

No Fan