

SHAKING WATER BATH

FEATUREA

- 1) The grounded circuit prevents a short circuit and overflow of electricity.
- 2) Equipped oveheat protection sensor.
- 3) Urethane foam optimizes insulation and heat loss
- 4) Time set fuction and constant temperature.
- 5) Uniquely designed exterior and mono piec bath.
- 6) Equipped internal bottom for the safety of the unit.
- 7) LED indicates tempreature, time and state of pwer statue
- 8) Digital timer fuctionW
- 9) Electronically controlled shaking mechanism minimizes the noise during working.



DISTINCTION

1. Separatable Spring Wire Multi Rack



Separatable Wire Rack maximizes customer's convenience for using it.

2. Stable Temp Control



The stable and accurate temp control is possible. by sensor.

3. Safety Switch



Safety switch control over temp for user and instrument.

SPECIFICATION

Spec. / Model		K-SWB25	K-SWB35	K-SWB55
Motion		Reciprocating		
Size	Internal (WxDxH, mm)	350x300x240	500x400x240	660x350x240
	External (WxDxH, mm)	750x365x430	720x365x430	880x420x430
Temp	Range	Ambient +5 °C ~ 95 °C		
	Controller	Digital PID Controller		
	Accurate	+/-0.2°C at 37°C		
Rpm	Range	20 to 180 rpm		
	Controller	Solid State Electronic Control, Digital Display		
	Stroke	25mm		
Timer		Digital Timer 99hr 59min or infinity		
Material	Internal	Stainless steel, Polishing Plate		
	External	Steel, Power Coating		
Capacity		25 liter	35 liter	55 liter
Wire		Spring wire Rack		

* Customized Model is possible.

** Above speficiation can be modified without notice.